



Product Datasheet

ERLIN1 antibody (orb382694)

Description	Rabbit polyclonal antibody to ERLIN1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human SPFH1. The exact sequence is proprietary.
Target	ERLIN1; ERLIN2
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	O94905 , O75477 , Q91X78 , B5DEH2 , Q8BFZ9
Entrez	244373 , 10613 , 290823 , 226144 , 11160
Dilution Range	WB: 1:500:2000, IHC-P: 1:50:200
Expiration Date	12 months from date of receipt.



Western blot analysis of SPFH1 expression...



Immunohistochemical analysis of SPFH1 staining...