


Product Datasheet

Collagen I antibody (orb382699)

Description	Rabbit polyclonal antibody to Collagen I
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human Collagen 1 alpha 1/2. The exact sequence is proprietary.
Target	COL1A1; COL1A2
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P02452 , P08123 , P02466 , Q01149 , P11087 , P02454
Entrez	1278 , 12842 , 12843 , 84352 , 1277 , 29393
Dilution Range	WB: 1:500:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of Collagen 1 alph...