

Product Datasheet

VEGF antibody (orb382731)

Description	VEGF antibody
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant protein of human VEGF
Target	VEGF
Form/Appearance	PBS with 0.05% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified by protein G.
Uniprot ID	P15692, Q00731, P16612
NCBI	001020537, 001025366
Entrez	5788
Dilution Range	WB: 1:100-500, IHC-P: 1:100-300
Expiration Date	12 months from date of receipt.