

Product Datasheet

CLU antibody (orb382739)

Description	CLU antibody
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant protein of human CLU
Target	CLU
Form/Appearance	PBS with 0.05% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified by protein G.
Uniprot ID	P05371, P01009, Q06890
NCBI	001822, 001831
Entrez	1191
Dilution Range	WB: 1:100-500, IHC-P: 1:100-300
Expiration Date	12 months from date of receipt.