

Product Datasheet

ICAM1 antibody (orb382743)

Description	ICAM1 antibody
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant protein of human ICAM1
Target	ICAM1
Form/Appearance	PBS with 0.05% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified by protein G.
Uniprot ID	P13597, P05362, Q00238
NCBI	000201, 000192
Entrez	3383
Dilution Range	WB: 1:100-500, IHC-P: 1:100-300
Expiration Date	12 months from date of receipt.