

Product Datasheet

FGFR1 antibody (orb38282)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	FGFR1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	91868
Uniprot ID	P11362
NCBI	NP_001167535.1 , NP_075598.2 , NP_075594.1 , NP_075593.1 , NP_001167536.1 , NP_056934.2 , NP_001167534.1 , NP_001167538.1 , NP_001167537.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

