




Product Datasheet

FGFR2 antibody (orb38286)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	FGFR2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	92025
Uniprot ID	P21802
NCBI	NP_001138390.1, NP_001138389.1, NP_001138391.1, NP_000132.3, NP_001138385.1, NP_075418.1, NP_075259.4, NP_001138386.1, NP_001138388.1, NP_001138387.1
Dilution Range	WB: 1:1000, IHC-P: 1:50-100, FACS: 1:10-50, IF/ICC: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of NCI-H460 cell l...



Immunofluorescence analysis of HeLa cell...