


Product Datasheet

Cystatin C protein (orb382860)

Description	Cystatin C protein
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, SDS-PAGE, WB
Target	Cystatin C
Form/Appearance	20mM Tris,150mM NaCl,0.5mM DTT 5% Glycerol,pH8.0
Concentration	0.5 mg/ml
Storage	Store at -20°C for at least 1 year. After reconstitution, protein solution is stable at -20°C for 6 months to 12 months. Keep at -80°C for long term storage. Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	N-Trx-His6
Purity	> 90%
MW	31.3 kDa
NCBI	000099
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Cystatin C protein