

Product Datasheet

UNC5B protein (orb382862)

Description	UNC5B protein
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, SDS-PAGE, WB
Target	UNC5B
Form/Appearance	50mM Tris,50mM NaCl,2M Urea,1mM DTT,pH8.0
Concentration	0.5 mg/ml
Storage	Store at -20°C for at least 1 year. After reconstitution, protein solution is stable at -20°C for 6 months to 12 months. Keep at -80°C for long term storage. Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	N-Trx-His6
Purity	> 80%
MW	42.4 kDa
NCBI	170744
Expiration Date	6 months from date of receipt.