


Product Datasheet

IL12B (phospho-40) protein (orb382869)

Description	IL12B (phospho-40) protein
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, SDS-PAGE, WB
Target	IL12B (phospho-40)
Form/Appearance	50mM Tris, 50mM NaCl, pH8.0
Concentration	0.5 mg/ml
Storage	Store at -20°C for at least 1 year. After reconstitution, protein solution is stable at -20°C for 6 months to 12 months. Keep at -80°C for long term storage. Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	N-His6
Purity	> 85%
MW	38.5 kDa
NCBI	067499
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of IL12B (phospho-40) ...