


Product Datasheet

MCP1 protein (orb382888)

Description	MCP1 protein
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, SDS-PAGE, WB
Target	MCP1
Form/Appearance	10mM PBS,pH7.4
Concentration	0.5 mg/ml
Storage	Store at -20°C for at least 1 year. After reconstitution, protein solution is stable at -20°C for 6 months to 12 months. Keep at -80°C for long term storage. Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	N-his6
Purity	> 90%
MW	12.5 kDa
NCBI	002982
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
CCL2 protein