


Product Datasheet

CXCL13 protein (orb382893)

Description	CXCL13 protein
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, SDS-PAGE, WB
Target	CXCL13
Form/Appearance	50mM Tris, 50mM NaCl, 5% Glycerol pH 8.0
Concentration	0.5 mg/ml
Storage	Store at -20°C for at least 1 year. After reconstitution, protein solution is stable at -20°C for 6 months to 12 months. Keep at -80°C for long term storage. Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	N-Trx-his6
Purity	> 90%
MW	28.3 kDa
NCBI	006419
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of CXCL13 protein