


## Product Datasheet

### IL6-M protein (orb382898)

|                            |                                                                                                                                                                                                  |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | IL6-M protein                                                                                                                                                                                    |
| <b>Reactivity</b>          | Mouse                                                                                                                                                                                            |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                                     |
| <b>Tested Applications</b> | ELISA, SDS-PAGE, WB                                                                                                                                                                              |
| <b>Target</b>              | IL6-M                                                                                                                                                                                            |
| <b>Form/Appearance</b>     | 50mM Tris, 50mM NaCl, pH8.0                                                                                                                                                                      |
| <b>Concentration</b>       | 0.2 mg/ml                                                                                                                                                                                        |
| <b>Storage</b>             | Store at -20°C for at least 1 year. After reconstitution, protein solution is stable at -20°C for 6 months to 12 months. Keep at -80°C for long term storage. Avoid repeated freeze/thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                                                                            |
| <b>Application notes</b>   | N-his6                                                                                                                                                                                           |
| <b>Purity</b>              | > 90%                                                                                                                                                                                            |
| <b>MW</b>                  | 25.5 kDa                                                                                                                                                                                         |
| <b>NCBI</b>                | <b>031168</b>                                                                                                                                                                                    |
| <b>Expiration Date</b>     | 6 months from date of receipt.                                                                                                                                                                   |



SDS-PAGE  
analysis of  
IL6-M protein