


Product Datasheet

IL27A-M protein (orb382900)

Description	IL27A-M protein
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, SDS-PAGE, WB
Target	IL27A-M
Form/Appearance	50mM Tris, 50mM NaCl, pH8.0
Concentration	0.5 mg/ml
Storage	Store at -20°C for at least 1 year. After reconstitution, protein solution is stable at -20°C for 6 months to 12 months. Keep at -80°C for long term storage. Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	N-his6
Purity	> 85%
MW	27.1 kDa
NCBI	145636
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of IL27A-M protein