


Product Datasheet

CXCL9-M protein (orb382901)

Description	CXCL9-M protein
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, SDS-PAGE, WB
Target	CXCL9-M
Form/Appearance	10mM PBS,pH7.4
Concentration	0.5 mg/ml
Storage	Store at -20°C for at least 1 year. After reconstitution, protein solution is stable at -20°C for 6 months to 12 months. Keep at -80°C for long term storage. Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	N-Trx-his6
Purity	> 90%
MW	30.1 kDa
NCBI	008599
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
CXCL9-M
protein