


Product Datasheet

IL4-M protein (orb382905)

Description	IL4-M protein
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, SDS-PAGE, WB
Target	IL4-M
Form/Appearance	10mM PBS,pH7.4
Concentration	0.5 mg/ml
Storage	Store at -20°C for at least 1 year. After reconstitution, protein solution is stable at -20°C for 6 months to 12 months. Keep at -80°C for long term storage. Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	N-Trx-his6
Purity	> 90%
MW	30.9 kDa
NCBI	021283
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
IL4-M protein