


Product Datasheet

14-3-3 zeta/delta antibody (orb382906)

Description

14-3-3 zeta/delta antibody


Species/Host

Rabbit

Reactivity

Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

ICC, IF, IHC-P, WB

Applications

Immunohistochemical staining of paraffin...

Immunogen

KLH conjugated synthetic peptide derived from human 14-3-3 zeta/delta. Please contact us for the exact immunogen sequence. The peptide is available as orb13048.

Target

14-3-3 zeta/delta

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Immunohistochemical staining of paraffin...

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

IHC-P staining of pig uterus tissue usin...

MW

28 kDa

Uniprot ID
[P63102](#), [P63101](#), [P63104](#)
NCBI
[001135699](#), [001129171](#)
Entrez
[7534](#)
Dilution Range

IF/ICC: 1:100-800, IHC-P: 1:100-800, WB: 1:200-2000

Expiration Date

12 months from date of receipt.

