

Product Datasheet

BNP antibody (orb382908)

Description

BNP antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ICC, IF, IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human BNP. Please contact us for the exact immunogen sequence. The peptide is available as orb13050.

Target

BNP

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Application notes

Not suitable for WB in Rat tissue

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

15 kDa

Uniprot ID

[P13205](#), [P40753](#), [P16860](#)

NCBI

[002512](#), [002521](#)

Entrez



[4879](#)

Dilution Range


IF/ICC: 1:100-800, IHC-P: 1:100-800, WB: 1:200-1500

Expiration Date


12 months from date of receipt.



IHC-P image of rat brain tissue using BNP...



Immunohistochemical staining of rat brain...



IF image of rat brain tissue using BNP a...