

Product Datasheet

Caspase 3 antibody (orb382909)

Description

Caspase 3 antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human Caspase 3. Please contact us for the exact immunogen sequence. The peptide is available as orb13051.

Target

Caspase 3

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

32 kDa

Uniprot ID

[P70677](#), [P42574](#), [P55213](#)

NCBI

[004337](#), [004346](#)

Entrez



[836](#)

Dilution Range


IHC-P: 1:30-300, WB: 1:200-1500

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of paraffin...



IHC-P staining of rat kidney tissue usin...



Western blot analysis of rat lymph node ...