

Product Datasheet

DSPP antibody (orb382912)

Description

DSPP antibody

Species/Host

Rabbit

Reactivity

Guinea pig, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human DSPP. Please contact us for the exact immunogen sequence. The peptide is available as orb13054.

Target

DSPP

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

131 kDa

Uniprot ID

[Q9NZW4](#), [Q62598](#), [P97399](#)

NCBI

[014208](#), [055023](#)

Entrez


[1834](#)

Dilution Range


IHC-P: 1:50-500, WB: 1:200-2000

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of paraffin...



Immunohistochemical staining of guinea p...



IHC-P staining of rat brain tissue using...