

Product Datasheet

EPX antibody (orb382913)

Description

EPX antibody

Species/Host

Rabbit

Reactivity

Guinea pig, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human EPX. Please contact us for the exact immunogen sequence. The peptide is available as orb13055.

Target

EPX

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

81 kDa

Uniprot ID

[D3ZSY4](#), [O88933](#), [P11678](#)

NCBI

[000493](#), [000502](#)

Entrez


[8288](#)

Dilution Range


IHC-P: 1:100-500, WB: 1:200-1000

Expiration Date


12 months from date of receipt.



Western blot analysis of rat lymph node ...



Immunohistochemical staining of paraffin...



IHC-P image of pig uterus tissue using E...