

Product Datasheet

MET antibody (orb382916)

Description

MET antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

ICC, IF, IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human MET. Please contact us for the exact immunogen sequence. The peptide is available as orb13058.

Target

MET

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

156 kDa

Uniprot ID

[P97523](#), [P16056](#), [P08581](#)

NCBI

[000236](#), [000245](#)

Entrez


[4233](#)

Dilution Range


IF/ICC: 1:100-500, IHC-P: 1:100-500, WB: 1:200-1500

Expiration Date


12 months from date of receipt.



Immunofluorescence analysis of rat pancreas tissue



IHC-P image of rat pancreas tissue using MET antibody



Immunohistochemical staining of paraffin-embedded rat pancreas tissue