

Product Datasheet

NKX6-1 antibody (orb382917)

Description

NKX6-1 antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human NKX6-1. Please contact us for the exact immunogen sequence. The peptide is available as orb13059.

Target

NKX6-1

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

38 kDa

Uniprot ID

[O35762](#), [P78426](#), [Q99MA9](#)

NCBI

[006159](#), [006168](#)

Entrez



[4825](#)

Dilution Range


IHC-P: 1:100-600, WB: 1:200-1500

Expiration Date


12 months from date of receipt.



IHC-P staining of rat stomach tissue using NKX6-1 antibody.



Western blot analysis of mouse stomach tissue using NKX6-1 antibody.



Immunohistochemical staining of paraffin-embedded tissue using NKX6-1 antibody.