



Product Datasheet

VAMP7 antibody (orb382918)


Description	VAMP7 antibody
Species/Host	Rabbit
Reactivity	Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human VAMP7. Please contact us for the exact immunogen sequence. The peptide is available as orb13060.
Target	VAMP7
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	25 kDa
Uniprot ID	Q9JHW5 , P70280 , P51809
Entrez	6845
Dilution Range	IHC-P: 1:100-500, WB: 1:200-1500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat stom...



IHC-P staining of rat spinal cord tissue...



Western blot analysis of rat lymph node ...