


Product Datasheet

SNX5 antibody (Biotin) (orb382924)

Description	Goat polyclonal antibody against SNX5 conjugated to Biotin.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	SNX5, Biotinylated
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 48kDa band observed in lysates of cell line A549 (calculated MW of 46.8kDa according to NP_055241.1). See non-biotinylated parental product's datasheet for further QC data. Recommended concentration: 1-3
Protein Sequence	SLLQSCIDLFKNN.
Clonality	Polyclonal
MW	46.8; 34.9
Entrez	27131
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot of A549 lysate using SNX5 a...