

Product Datasheet

KPNA3 antibody (Biotin) (orb382925)

Description Goat polyclonal antibody against KPNA3 conjugated to Biotin.

Species/Host Goat

Reactivity Bovine, Canine, Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ELISA, ICC, IF, IHC, IP, WB

Target KPNA3 / IPOA4, Biotinylated

Preservatives Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage Aliquot and store at -20°C. Minimize freezing and thawing.

Note For research use only

Application notes WB: Approx 55kDa band observed in Rat Testis lysates (calculated MW of 57.8kDa according to Human NP_002258.2 and Rat NP_001014792.1). See non-biotinylated parental product's datasheet for further QC data. Recommended concentration: 1-3

Protein Sequence CDPTANLQTKEFNF.


Clonality Polyclonal

MW 57.8

Entrez [3839](#)

Dilution Range ELISA: 1:32000, WB: 1-3 µg/ml

Expiration Date 12 months from date of receipt.



Western blot of Rat Testis lysate using

...