


Product Datasheet

PAX3 antibody (Biotin) (orb382926)

Description	Goat polyclonal antibody against PAX3 conjugated to Biotin.
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	PAX3, Biotinylated
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 45kDa band observed in Human Duodenum lysates (calculated MW of 45.2kDa according to NP_852125.1 and 44.8kDa according to NP_852126.1). See non-biotinylated parental product's datasheet for further QC data. Recommended concentration: 0.3-1
Protein Sequence	TTLAGAVPRMMC.
Clonality	Polyclonal
MW	53.0; 24.1; 22.7; 53.5; 56.0; 44.8; 45.2; 53.3
Entrez	5077
Dilution Range	ELISA: 1:32000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot
of Human
Duodenum
lysate us...