


Product Datasheet

Nsg2 antibody (orb382931)

Description	Goat polyclonal antibody to Nsg2
Species/Host	Goat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Target	Nsg2 (mouse, C terminus)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 19kDa band observed in Mouse fetal Brain lysates (calculated MW of 19.0kDa according to NP_032767.1). Recommended concentration: 0.1-0.3
Protein Sequence	HEPKPPKTQGH
Clonality	Polyclonal
MW	19.0; 11.3
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse fetal Bra...