

Product Datasheet

GSTA4 antibody (orb382934)

Description	Goat polyclonal antibody to GSTA4
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	GSTA4
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 25kDa band observed in Human Liver lysates (calculated MW of 25.7kDa according to NP_001503.1). Recommended concentration: 0.3-1
Protein Sequence	QEYTVKLSNIPTIK
Clonality	Polyclonal
MW	25.7
Entrez	2941
Dilution Range	ELISA: 1:64000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Liver lys...