

Product Datasheet

GFAP antibody (orb382936)

Description	Goat polyclonal antibody to GFAP. Glial fibrillary acidic protein is a major intermediate filament p...
Species/Host	Goat
Reactivity	Bovine, Canine, Donkey, Equine, Feline, Gallus, Goat, Guinea pig, Hamster, Human, Mouse, Other, Porcine, Primate, Rabbit, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant peptide derived from within residues 375 aa to the C-terminus of human GFAP produced in E. coli.
Target	GFAP
Preservatives	PBS, 20% glycerol and 0.05% sodium azide
Form/Appearance	Polyclonal antibody supplied as a 100 µl (3 mg/ml) aliquot in PBS, 20% glycerol and 0.05% sodium azide. This antibody is epitope-affinity purified from goat antiserum.
Concentration	3 mg/ml
Storage	Store at -20°C for long-term storage. Store at 2-8°C for up to one month.
Note	For research use only
Application notes	The antibody solution should be gently mixed before use.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB:1:500-1:5,000
Expiration Date	12 months from date of receipt.