

Product Datasheet

Human DLAT DLTA protein (orb382949)

Description	Recombinant human DLAT DLTA protein.
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	E.coli and Yeast N-terminal 6xHis-tagged Partial
Protein Sequence	KVPLPSLSPTMQAGTIARWEKKEGDKINEGDLIAEVETDKATVGFESLEECYMAKI LVAEGTRDVPIGAIICITVGKPEDIEAFKNYTLDSSAATPQAAPAPTPAATASPPTP SAQAPGSSYPHMQVLLPALSPTMTMGTVQRWEKKVGEKLSEGDLLAEIETDKA TIGFEVQEEGYLAKILVPEGTRDVPLGTPLCIIVEKEADISAFADYRPTEVTDLKPQ VPPPTPPPVAAPPTPQPLAPTPSAPCPATPAGPKGRVFSPLAKKLAVEKGIDLT QVKGTGPDGRITKKDIDSFVPSKVAPAPAAVPPPTGPGMAPVPTGVFTDIPISNIR RVIAQRLMQSKQTIPHYLSIDVNMGEVLLVRKELNKILEGRSKISVNDFIKASAL ACLKVPEANSSWMDTVIRQNHVVDVSVAVSTPAGLITPIVFNAHIKGVETIANDV VSLATKAREGKLQPHEFQGGTFTISNLGMFGIKNFSAINPPQACILAIGASEDKLV PADNEKGFDAVSMMSVTLSCDHRVVDGAVGAQWLAEFRKYLE
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	62.2 kDa
Uniprot ID	P10515
Expiration Date	6 months from date of receipt.