

Product Datasheet

Human NKRF protein (orb382964)

Description

Recombinant human NKRF protein.

Endotoxins

Not test.

Preservatives

If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.

Form/Appearance

Liquid or Lyophilized powder

Storage

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Tag

N-terminal 6xHis-tagged

Note

For research use only

Application notes

E.coli and Yeast N-terminal 6xHis-tagged Partial

Protein Sequence

ILQMAEGIDIGEMPSYDLVLSKPSKGQKRHLSTCDGQNPPKKQAGSKFHARPRFE
PVHFVASSSKDERQEDPYGPQTKEVNEQTHFASMPRDIYQDYTDQSFSDGNS
QYCDSSGFILTKDQPV TANMYFDSGNPAPSTTSQQANSQSTPEPSQTFPESV
VAEKQYFIEKLTATIWKNSNPEMTSGSDKINYTYMLTRCIQACKTNPEIYAPLKE
IPPADIPKNNKLLTDGYACEVRCQNIYLTGTYAGSKNGSRDRATELAVKLLQKRIE
VRVVRKFKHTFGEDLVVCQIGMSSYEFPPALKPPEDLVVLGKDASGQPIFNASA
KHWTFNVITENANDAIGILNNSASFNMKSIEYKEMMPNRTWRCRVFLQDHCLA
EGYGTKKTSKHAAADEALKILQKTQPTYPSVKSSQCHTGSSPRGSGKKKDIKDLV
VYENSSNPVCTLNDTAQFNRMTVEYVYERMTGLRWKCKVILESEVIAEAVGVKK
TVKYEAAGEAVKTLKKTQPTVINNLKKGAVEDVISRNEIQGRSAEEAYKQKIKED
NIGNQLLRKMGWTGGGLGKSGEGIREPISVKEQHKREGLDVERVNIKARDI
EQIIRNYARSESHDTLTF SRELTNDERKQIHQIAQKYGLKSKSHGVGHDRYLVVG
RKRRKEDLLDQLKQEGQVGHYELVMPQ

Purity

Greater than 90% as determined by SDS-PAGE.

Source

E.coli

MW

81.1 kDa

Uniprot ID
015226
Expiration Date

6 months from date of receipt.