




Product Datasheet

VEGFR2 antibody (orb38298)

Description	Rabbit polyclonal antibody to VEGFR2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	KDR
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	151527
Uniprot ID	P35968
NCBI	NP_002244.1
Dilution Range	WB: 1:1000, IF/ICC: 1:10-50, FACS: 1:10-50, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Immunohistochemical staining of paraffin...