

Product Datasheet

Human ESR1 protein (orb382994)

Description	Recombinant human ESR1 protein.
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	E.coli and Yeast N-terminal 6xHis-tagged Partial
Protein Sequence	SGMALLHQIQGNELEPLNRPQLKIPLERPLGEVYLDSSKPAVYNYPEGAAYEFNA AAAANAQVYGQTGLPYGPGSEAAAFGSNGLGGFPPLNSVSPSPLMLLHPPQLS PFLQPHGQQVPYYLENESGYTVREAGPPAFYRPNSDNRRQGGRERLASTNDK GSMAMESAKETRYCAVCNDYASGYHYGVWSCEGCKAFFKRSIQGHNDYMCPA TNQCTIDKNRRKSCQACRLRKCYEVGMMKGGIRKDRRGGRMLKHKRQRDDGE GRGEVGSAGDMRAANLWPSPLMIKRSKNSLALSLTADQMVSALLDAEPPILYS EYDPTPRPFSEASMMGLLTNLADRELVHMINWAKRVPGFVDLTLHDQVHLLECA WLEILMIGLVWRSMHPGKLLFAPNLLLDNRNQGKCEGMVEIFDMLLATSSRFR MMNLQGEFVCLKSIILLNSGVYTFLSSTLKSLEEKDHIHRVLDKITDTLIHLMKA GLTLQQQHQRQAQLLLILSHIRHMSNKGMEHLYSMKCKNVVPLYDLLLEMLDAH RLHAPTSRGGASVEETDQSHLATAGSTSSHSLQKYITGEAEGFPATV
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Yeast
MW	67.2 kDa
Uniprot ID	P03372
Expiration Date	6 months from date of receipt.