

Product Datasheet

Human OGN protein (orb383003)


Description	Recombinant human OGN protein.
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid: Tris/PBS-based buffer, 5%-50% glycerol.
Storage	Store at -20°C/-80°C upon receipt, aliquoting is necessary for mutiple use. Avoid repeated freeze-thaw cycles.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	E.coli and Yeast N-terminal 6xHis-tagged Full Length
Protein Sequence	PPTQQDSRIIYDYGTDNFEESIFSQDYEDKYLDGKNIKEKETVIIPNEKSLQLQKDE AITPLPPKKENDEMPPTCLLCVCLSGSVYCEEVDIDAVPPLPKE SAYLYARFNKIKKLTAKDFADIPNLRRLDFTGNLIEDIEDGTFSKLSLLEELSLAEN QLLKLPVLPPKLTTFNAKYNKIKSRGIKANAFKLLNLTFLYLDH NALESVPLNLPESLRVIHLQFNNIASITDDTFCKANDTSYIRDRIEEIRLEGNPVLG KHPNSFICKRLPIGSYF
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Yeast
MW	33.7 kDa
Uniprot ID	P20774
Expiration Date	6 months from date of receipt.



Based on the
SEQUEST from
database of
Ye...



Based on the
SEQUEST from
database of
Ye...



SDS-PAGE
(reduced) with
5% enrichment
ge...