

Product Datasheet

Hamster P4HB protein (orb383096)

Description	Recombinant hamster P4HB protein.
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	E.coli and Yeast N-terminal 6xHis-tagged Full Length
Protein Sequence	DAPEEEDNVLVLKSNFAEALAAHNYLLVEFYAPWCGHCKALAPEYAKAAAKLK AEGSEIRLAKVDATEESDLAQYGVGYPTIKFFKNGDTASPKEYTAGREADDIV NWLKKRTGPAATTLSDTAAAETLIDSSEVAVIGFFKDVESDSAKQFLAAEAVDDI PFGITSNSGVFSKYQLDKDGVVLFKKFDEGRNNFEGEVTKELLDIFIKHNQLPLVI EFTEQTAPKIFGGEIKTHILLFLPKSVSDYDGKLGNFKAAEGFKGKILFIFIDSDHT DNQRILEFFGLKKEECPAVRLITLEEEMTKYPESDELTAEKITEFCHRFLEGKIKP HLMSQELPEDWDKQPVKVLVGKNFEEVAFDEKKNVFEFYAPWCGHCKQLAPI WDKLGETYKDHEIIIAKMDSTANEVEAVKVSFPTLKFFPATADRTVIDYNGER TLDGFKKFLESQGDGAGDDDDVDLEEALPDMEEDDDQKAVKDEL
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	59 kDa
Uniprot ID	Q8R4U2
Expiration Date	6 months from date of receipt.