


Product Datasheet

Amphibian catesbeiana Saxiphilin protein (orb383116)

Description	Recombinant amphibian catesbeiana Saxiphilin protein.
Reactivity	Amphibian
Target	N/A
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	• E.coli and Yeast N-terminal 6xHis-tagged • Partial
Protein Sequence	AHLPSKNKVRWCTINKLEKMKCDDWSAVSGGAIACEASCPKGCVKQILKGEAD AVKLEVQYMYEALMCGLLPAVEEYHNKDDFGPCKTPGSPYTDFTGLRAVALVKK SNKDINWNNIKGKKSCHTGVGDIAGWVIPVSLIRRQNDNSDIDSFFGESCAPGS DTKSNLCKLCIGDPKNSAANTKCSLSDKEAYYGNQGAFRCLVEKGDVAFVPHTV VFENTDGKNPAVWAKNLKSEDFELLCLDGSRAPVSNYKSCKLSGIPPPAIVTREE SISDVVRIVANQQSLYGRKGFEKDMFQLFSSNKGNNLLFNDNTQCLITFDRQPK DIMEDYFGKPYTTTVYGASRSAMSELISACTIKHC
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Yeast
MW	41.7 kDa
Uniprot ID	P31226
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
using
Amphibian
cat...