

Product Datasheet

Human SIRPG protein (orb383170)

Description	Recombinant human SIRPG protein.
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	E.coli and Yeast N-terminal 6xHis-tagged Extracellular Domain
Protein Sequence	EEELQMIQPEKLLLVTVGKTATLHCTVTSLLPVGPVLWFRGVGPGRELIYNQKEG HFPRVTTVSDLTKRNNMDFSISSITPADVGTYCVKFRKGSPENVEFKSGPGTE MALGAKPSAPVVLGPAARTTPEHTVSFTCESHGFSPRDITLKWFKNGNELSDFQ TNVDPTGQSVAYSIRSTARVVLDPWDVRSQVICEVAHVTLQGDPLRGTTANLSEA IRVPPTLEVTQQPMRVGNQVNVTCQVRKFYPQSLQLTWSSENGNVCQRETASTL TENKDGTYNWTWFLVNIQDQRRDDVLTQVKHGDGLAVSKRLALEVTVHQQ DQSSDATP
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Yeast
MW	38.8 kDa
Uniprot ID	Q9P1W8
Expiration Date	6 months from date of receipt.