

Product Datasheet

ILK antibody (orb38319)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ILK
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	51419
Uniprot ID	Q13418
NCBI	NP_004508.1 , NP_001014794.1 , NP_001265370.1 , NP_001014795.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

