

Product Datasheet

E. coli dnaK protein (orb383248)

Description	Recombinant e. coli dnaK protein.
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	E.coli and Yeast N-terminal 6xHis-SUMO-tagged Full Length
Protein Sequence	MGKIIIGIDLGTNSCVAIMDGTTPRVLENAEGDRTTPSIIAYTQDGETLVGQPAKR QAVTNPQNTLFAIKRLIGRRFQDEEVQRDVSIMPFKIIAADNGDAWVEVKGQKM APPQISAEVLKKMKKTAEDYLGEPVTEAVITVPAYFNDAQRQATKDAGRIAGLEV KRIINEPTAAALAYGLDKGTGNRTIAVYDLGGGTFDISIIEIDEVDGEKTFEVLATN GDTHLGGEDFDSRLINYLVEEFKKDQGIDLRNDPLAMQRLKEAAEKAKIELSSAQ QTDVNLPIYTADATGPKHNMNIKVTRAKLESLVEDLVNRSIEPLKVALQDAGLSVS DIDDVILVGGQTRMPMVQKKVAEFFGKEPRKDVNPDEAVAIGA AVQGGVLTGD VKDVLLLDVTPLSLGIETMGGVMTTLIAKNNTIPTKHSQVFSTAEDNQSAVTIHVL QGERKRAADNKS LGQFNLDGINPAPRGMPQIEVTFDIDADGILHVS AKDKNSGK EQKITIKASSGLNEDEIQKMVRDAEANA EADRKFEELVQTRNQGDHLLHSTRKQ VEEAGDKLPADDKTAIESALTALETALKGEDKAAIEAKMQELA QVSQKLMEIAQQ QHAQQQTAGADASANNAKDDDDVDAEFEEVKDKK
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	73.1 kDa
Uniprot ID	Q1RGI8
Expiration Date	6 months from date of receipt.