




Product Datasheet

c-Kit antibody (orb38330)

Description	Rabbit polyclonal antibody against c-Kit
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	KIT
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	109865
Uniprot ID	P10721
NCBI	NP_000213.1 , NP_001087241.1
Dilution Range	WB: 1:1000, FACS: 1:10-50, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa whole cell...



Flow cytometric analysis of NCI-H460 cel...