

Product Datasheet

Pig GPX5 protein (orb383352)

Description	Recombinant pig GPX5 protein.
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged and C-terminal Myc-tagged
Note	For research use only
Application notes	E.coli and Yeast N-terminal 6xHis-tagged Full Length of Mature Protein
Protein Sequence	NSNLEKMDCYKDVTGTIYDYDAFTLNGNEHIQFKQYAGKHLFVNVATYCGLTA QYPELNTLQEELKPFGLVVLGFPCNQFGKQEPGENSEILLGLKYVRPGGGYVPNF QLFEKGDVNGEKEQKVFTFLKHSCPHPSSELIGSIGYISWEPIRVHDIRWNFEKFLV GPDGVPVMRWVHETPISTVKSDILAYLKQFKTE
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Yeast
MW	26.1 kDa
Uniprot ID	O18994
Expiration Date	6 months from date of receipt.