

Product Datasheet

Birds SP protein (orb383393)

Description	Recombinant birds SP protein.
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	Tag-Free
Note	For research use only
Application notes	E.coli and Yeast NO-tagged Full Length of Mature Protein
Protein Sequence	WEWPWNRKPTKFIPISPNPRDKWCRLNLGPAWGGRC
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Yeast
MW	4.3 kDa
Uniprot ID	P05623
Expiration Date	6 months from date of receipt.