

Product Datasheet

Bacterial egsA protein (orb383426)

Description	Recombinant bacterial egsA protein.
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	E.coli and Yeast N-terminal 6xHis-tagged Full Length
Protein Sequence	MNQTL EEV VSLARQC GCGHRHYDIP IEQM VVGREAFARLVAYLRHKRYERVAIV ADDHTFAAVERS LCDQLENGSIRYTVCLVQPDENG DVIADERSIVQV LLETPDDV DVLIAGAGTIHDITRFSSYKM RIPFISVPTAPSV DGTSMGAPLIIRGVK KTIQAQA PIAVFAHTGVLCQSPKEMIAAGFGDMVAKYTS LADWQFAHWMAD EPCPFVHQ LTEQSLQTCVDHIDDIAAGGEQGIRVLMDALLQSGIAMLLMGQSY SASGAEHHL SHYWEMEFLRQKKRQVLHGAKVGVSTPIII EHYQRVFWPLLNELEKRPKSMDEAT WERLKANTASIRELLES LPSPERIRTMLAKVGGAI APEQLGIDPQLVERSLRE AHR LRLNRFTMLYFLNELIFVE
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Yeast
MW	47.4 kDa
Uniprot ID	B3EYN8
Expiration Date	6 months from date of receipt.