

Product Datasheet

Mouse PCMP-H63 protein (orb383430)

Description	Recombinant mouse PCMP-H63 protein.
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Note	For research use only
Application notes	E.coli and Yeast N-terminal 10xHis-tagged and C-terminal Myc-tagged Full Length of Mature Protein
Protein Sequence	QNWQSSDGCSSYGTTGNGVPQENNTGGNHFQQDHS GHSSLD ELDSICREGKV KKAVEIIKSWRNEGYVVDLPRLFWIAQLCGDAQALQEAKVVHEFITSSVGISDISA YNSIIMYSGCGSVEDALT VFN S M PERNLETWCGVIRCF AKNGQGEDAIDTFSRF KQEGNKPDGEMFKEIFFACGVLGDMNEGLLHFESMYKEYGIIPCMEHYVSLVKM LAEPGYLDEALRFVESMEPNVDLWETLMNLSRVHGD LILGDR CQDMVEQLDAS RLNKESKAGLVPVKSSDLVKEKLQRMAGPNYGIRYMAAGDISRPENRELYMAL KSLKEHMI EIGYVPLSKLALHDVDQESKDENLFNHNERFAFISTFLDTPARSLIRV MKNLRVCADCHNALKLMSKIVGRELISRD AKRFHHMKDGVCS CREYW
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	53.3 kDa
Uniprot ID	Q9SUU7
Expiration Date	6 months from date of receipt.