

Product Datasheet

Bacterial rgpA protein (orb383435)

Description	Recombinant bacterial rgpA protein
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	Baculovirus and Mammalian Cell N-terminal 6xHis-tagged Partial
Protein Sequence	YTPVEEKQNGRMIVIVAKKYEGLDFVDWKNQRGLRTEVKVAEDIASPVTANAI QQFVKQEYEKEGNDLTYVLLVGDHDKDIPAKITPGIKSDQVYGQIVGNDHYNEVFI GRFSCESKEDLKTQIDRTIHYERNITTEDKWLGGALCIASAEAGGPSADNGESDIQ HENVIANLLTQYGYTKIIKCYDPGVTPKNIIDAFNGGISLVNYTGHGSETAWGTSH FGTTHVKQLTNSNQLPFIFDVACVNGDFLFSMPCFEALMRAQKDGKPTGTVAII ASTINQSWASPMRGQDEMNEILCEKHPNNIKRTFGGVTMNGMFAMVEKYKKDG EKMLDTWTVFGDPSLLVRTLVP TKMQVTAPAQINLT DASVNVSCDYNGAIATISA NGKMFGSAVVENG TATINLTGLTNESTLTLTVVGYNKETVIKTINTNGEPNPYQPV SNLTATTGGQKVT LKWDAPSTKT NATTN TARSVDGIRELVLLSVSDAPELLRS
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Mammalian cell
MW	58 kDa
Uniprot ID	P28784
Expiration Date	6 months from date of receipt.