

Product Datasheet

Human CRBN protein (orb383449)

Description	Recombinant human CRBN protein
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	C-terminal 9xHis-tagged
Note	For research use only
Application notes	Baculovirus and Mammalian Cell N-terminal 6xHis-tagged Full Length
Protein Sequence	MAGEGDQQDAAHNMGNHLLPAESEEDEMEVEDQDSKEAKKPNINFDTSLS PTSHTYLGADMEEFHGRTLHDDDDSCQVIPVLPQVMMILIPGQTLPLQLFHPQEV MVRNLIQKDRTFAYLAYSINVQEREAQFGTTAEIYAYREEQDFGIEIVKVKAIQRQR FKVLELRTQSDGIQAKVQILPECVLPSTMSAVQLESLNKCKIFPSKPVSRDQCS YKWWQKYQKRKFHCANLTSWPRWLYSLYDAETLMDRIKKQLREWDENLKDD LPSNPIDFSYRVAACLPIDDLRIQLLKIGSAIQRLRCELDINMKCTSLCCKQCQET EITTKNEIFSLSLCGPMAAYVNPBGYVHETLTVYKACNLNLIGRPSTEHSWFPGYA WTVACQCKICASHIGWKFTATKKDMSPQKFWGLTRALLPTIPDTEDEISPDKVILC L
Purity	Greater than 85% as determined by SDS-PAGE.
Source	Baculovirus
MW	52.5 kDa
Uniprot ID	Q96SW2
Expiration Date	6 months from date of receipt.