

Product Datasheet

Pig HSP90AA1 protein (orb383455)

Description	Recombinant pig HSP90AA1 protein
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	Baculovirus and Mammalian Cell N-terminal 6xHis-tagged Partial
Protein Sequence	VEKERDKEVSDDEAEKEDKEEKEKEEKESEDKPEIEDVGSDEEEEEKKDGDGK KKKKKKIKEKYIDQEELNKTPIWTRNPDDITNEEYGEFYKSLTNDWEDHLAVKHFS VEGQLEFRALLFVPRRAPFDLFENRKKKNNIKLYVR
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Baculovirus
MW	19.7 kDa
Uniprot ID	O02705
Expiration Date	6 months from date of receipt.