

Product Datasheet

E. coli insE1 protein (orb383462)

Description	Recombinant e. coli insE1 protein
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	C-terminal hFc-tagged
Note	For research use only
Application notes	Baculovirus and Mammalian Cell N-terminal FC-tagged Full Length
Protein Sequence	MTKTVSTSKKPRKQHSPEFRSEALKLAERIGVTAAARELSLYESQLYNWRSKQQN QQTSSERELEMSTEIARLKRQLAERDEELAILQKAATYFAKRLK
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Mammalian cell
MW	37.5 kDa
Uniprot ID	P0CF66
Expiration Date	6 months from date of receipt.