

Product Datasheet

Mouse Il7r protein (orb383464)

Description	Recombinant mouse Il7r protein
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	C-terminal 9xHis-tagged
Note	For research use only
Application notes	Baculovirus and Mammalian Cell N-terminal 6xHis-tagged Partial
Protein Sequence	ESGNAQDGDLEDADADDHSFWCHSQLEVDGSQHLLTCAFNDSDINTANLEFQI CGALLRVKCLTLNKLQDIYFIKTSEFLIGSSNICVKLGQKNLTCKNMAINTIVKAEA PSDLKVYRKEANDFLVTFNAPHLKKKYLKKVKHDVAYRPARGESNWTHVSLFH TRTTIPQRKLRPKAMYIEIKVRSIPHNDYFKGFWSEWSPSSTFETPEPKNQGGWD
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Baculovirus
MW	27 kDa
Uniprot ID	P16872
Expiration Date	6 months from date of receipt.